

Making Portable Oxygen *easy*

Model	B 160 Litre Oxygen Cylinder	Rove 6 (with double battery)	SimplyGo Mini (with double battery)	CH 460 Litre Oxygen Cylinder	Zen-O (with single battery)	SimplyGo (NLA) (with single battery)
Weight kg	1.3 kg	3.3 kg	3.7 kg	5.5 kg	6 kg	8.2 kg
Duration Setting 2 Pulse / OCD	5 hours and 48 minutes	10 hours	9 hours	17 hours	4 hours	3 hours and 30 minutes
Duration Setting 2 Continuous	1 hour and 18 minutes	n/a	n/a	3 hours and 54 minutes	45 minutes	42 minutes
Duration Setting 6 Pulse / OCD	2 hours	2 hours and 18 minutes	Setting 5 = 4 hours	5 hours and 42 minutes	1 hour and 45 minutes	1 hour and 54 minutes

Calculations based on 20 breaths per minute

This is a guide only for prescribers; please take into consideration different patient needs. Changing from one mode to another may not be suitable in all situations.

Always read the instructions for use before operating a medical device.

Please do not hesitate to contact us for any assistance.

We will help you find the POC that meets your oxygen needs and complements your way of life.

How do I contact Air Liquide Healthcare?

 alhenquires@airliquide.com

 1300 36 02 02

 au.healthcare.airliquide.com/oxygen

 store.airliquidehealthcare.com.au